



SPG MEDISAFE PRIVATE LIMITED
Imagination | Innovation | Research



UROSAFE™

PRODUCT
CATALOGUE



Efficient • Transparent • Inclusive

URETERAL INDWELLING STENT

DJ STENT - TECO FLEX



Applications

Used for temporary internal drainage from the Ureteropelvic junction to the bladder.

Type: Normal DJ Stent & LongDuro DJ Stent

(*Caution: Periodic evaluation is advised.
Stent must not remain indwelling more than three months).

OEO - One end open

BEO - Both end open

Highly Radiopaque for better visualization

Perfect shape to prevent migration of stent

Tapered end for easy insertion over guidewire

Sizes

Size (in Fr):

2.5 / 3 / 3.5 / 3.8 / 4 / 4.5 / 4.8 / 5 / 5.5 / 6 / 6.5 / 7 / 8

Length (in cms):

8cm - 30cm

Recommended Guide Wire (in Inch) 0.018" to 0.038"

Other sizes can be provided on request

Also Available in Multiloop

ULTRA LONG DURO



Key Features

- › High Quality Implant Grade blended polymer Material
- › Highly Radiopaque for better visualization & very good ability
- › Can be indwelled in body up to 6 months

Sizes

Size (in Fr):

4.7 / 5 / 6

Length (in cms):

24cm - 26cm

ENDOPYELOTOMY STENT



Applications

Used for temporary internal drainage from kidney to urinary bladder.
Provides optimal healing post Endopyelotomy. Soft and flexible for patient comfort.

Sizes

Size (in Fr):	4/8	5/10	6/12	7/14
Length (in cms):	16-28			

CATHETER

URETERAL CATHETER



Applications

Ureteral Catheters are used to insert into the ureter to prevent obstruction of the urine flow. Catheters are made and design using finest quality polyurethane material.
Catheter inserted into the ureter, either through urethra and bladder or posteriorly via kidney, used for drainage, retrograde pyelogram, atraumatic of ureteral orifice & navigation of torturous ureter, ureteral dilatation.

Key Features

- ▶ Open end tip allows placement over guide-wire for negotiation of a looped or stenotic (torturous) ureter.
- ▶ Bulb tip configuration concludes orifice and stabilize catheter by temporarily wedging in.
- ▶ Angled tip allows the passage of the catheter pass a ureteral stone in preparation for ESWL smooth surface for easy insertion.

Variants

OE - One end	CE - Closed end	BT - Bulb Tip
--------------	-----------------	---------------

Sizes

Size (in Fr):	3FR / 4FR / 5FR / 6FR / 7FR
Length (in cms):	70

Recommended Guide Wire (in Inch) 0.018" to 0.038"
Other sizes can be provided on request

SUPRAPUBIC CATHETER SET



Applications

Used for percutaneous suprapubic urinary bladder catheterization to relieve urinary bladder retention.

Key Features:

- › Ideal for urinary bladder drainage is required in emergency situation.
- › Helpful when urethral catheterization is not possible.
- › Comprehensive procedural kit for easy use.

*Also available in complete set

Sizes

Size (in Fr):	8	10	12	14	16	18	20
Length (in cms):	25cm (Malecot) / 40cm (Foley)						

MALECOT CATHETER



Applications

Used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage. This is a winged type catheter like a mushroom with slots in it. These catheters are best for drainage of thick viscous fluids. It is generally used to provide drainage following open renal or bladder surgery.

Key Features

- › Highly Radiopaque for better visualization.
- › Malecot flower design enhance drainage and promote catheter retention.
- › Round tip malecot catheters reduces the trauma to the tissues.
- › Malecot flower design is ideally suited for the drainage of thick viscous fluid.

Sizes

Size (in Fr):	8 / 10 / 12 / 14 / 16 / 18 / 20 / 22 / 24
Length (in cms):	22cm / 30cm

Other sizes can be provided on request

Recommended Guide Wire (in Inch) 0.038" St tip
Percutaneous Malecot Nephrostomy Catheter Set

*Malecot Catheter + Initial Puncture Needle + Guide wire + Urine bag Connector + Fascial Dilators + Scalpel

PIGTAIL CATHETER



Applications

Used for the temporary or permanent drainage & irrigation of urine from kidney direct puncture in presence of distal obstruction..
Used for the temporary drainage of urine from kidney direct by percutaneous placement, in the presence of distal obstruction puncture in Pigtail Nephrostomy Catheter.

Key Features

- › Highly radiopaque for better visualization
- › Pigtail tip configuration for the retention of the distal tip in the renal pelvis.
- › Locking types has more retention and can be kept for no longer duration.
- › Soft Flexible material for patient's comfort.

Sizes

Size (in Fr):	6 / 8 / 10 / 12 / 14 / 16 / 18 / 20 / 22 / 24
Length (in cms):	22cm / 30cm

Other sizes can be provided on request

Recommended Guide Wire (in Inch) 0.038" J tip
Percutaneous Pigtail Nephrostomy Catheter Set

* Pigtail Catheter + IP Needle + Guide wire + Urine bag Connector + Facial Dilators + Scalpel

URETERAL ACCESS SHEATH CATHETER



Applications

The ureteral access sheath facilitates the use of flexible ureteroscopy, enabling improved minimally invasive management of complex upper urinary tract diseases.


Key Features:

- › Kink-resistant sheath is reinforced with a stainless steel coil, designed to provide stability and strength through tortuous ureters
- › Hydrophilic coating on external sheath and dilator tip designed for ease of insertion

Sizes

Size (in Fr):	9.5FR / 11.5 FR	10FR / 12FR
Length (in cms):	35 cm / 45 cm	

CLEAN INTERMITTENT CATHETER




Applications

Clean Intermittent Catheterisation (CIC): Refers to the act of draining urine by passing a catheter through the urethra, past the sphincter into the bladder. The catheter is removed after the urine has been drained. This process is performed using a clean technique.

Sizes

Size (in Fr):	8FR - 26FR
Length (in cms):	30 cm / 40 cm

NEPHROSTOMY TUBE



Applications

A nephrostomy is a tube that lets urine drain from the kidney through an opening in the skin on the back. A thin, flexible tube goes through the opening and into the kidney. This is called a nephrostomy tube or catheter. Urine drains through the tube into a bag outside your body.

Sizes

Size (in Fr):	8FR - 26FR
Length (in cms):	30 cm / 40 cm

DILATOR

URETERAL BALLOON DILATOR



Applications

The Ureteral balloon dilator is inflated slowly to allow for accommodation of the ureter and no more than 2 atmospheres per minute is used. After dilation the balloon is deflated completely and removed. An occlusion tip retrograde contrast medium study is done at the completion of each procedure to evaluate ureteral injury.

Sizes

Size (in Fr):	5mm x 40mm
Length (in cms):	75 cm

URETERAL DILATOR SET

Used for dilation of the Ureter prior to Ureteroscopy procedures.

Sizes

In Fr:	Length (in cms):
6 - 16	70

*Used with 0.038 Guide wire




URETHRAL DILATOR SET

Used for dilation of the Urethra prior to cystoscopy procedures.

Sizes

In Fr:	Length (in cms):
8 - 22	40

*Customizable
*Used with 0.038 Guide wire




FASCIAL DILATOR

Used for fascial dilation.

Sizes

In Fr:	Length (in cms):
6 - 16	22

*Customizable
*Used with 0.038 Guide wire



MEATAL DILATOR

Used for self dilation of urethral meatus

Sizes

In Fr:	Length (in cms):	Gender
8 - 30	9	Male
8 - 30	9	Female
6 - 16	9	Pediatric



SCREW DILATOR

Used for one step Dilation This can be used instead of sequential dilation.

Sizes

In Fr:	Length (in cms):
6 - 14	22
6 - 16	22

*Customizable

*Used with 0.038 Guide wire



AMPLATZ DILATOR SET

Amplatz / Renal dilations and sheath are used to dilate nephrostomy tracts to allow use of instruments for procedures.

Sizes

In Fr:	Length (in cms):
6 - 30	30



GUIDE WIRES

Guide wires are intended to facilitate the placement of endo urological Devices during diagnostic or interventional procedures. Used for gaining access and navigating torturous urinary tract.



SS GUIDE WIRES - PTFE COATED

Sizes

Size Length	"0.018 to 0.038" - 150cm
Type	'J' / 'St.' Tip

SLIDE GUIDE WIRES - NITINOL

Sizes

Size Diameter	"0.025 to 0.038" - 150cm
Type	'St' Tip

STRIPED GUIDE WIRES - ZEBRA


Sizes

Size Diameter	"0.018 to 0.038" - 150cm
Type	'St' Tip

*Also available in 260cm length

BASKET

LOOP GAGE BASKET



Sizes

In Fr: Length (in cms):

1.9 / 2.2

90 / 120

Accept 280 / 365 micron fiber

- › Nitinol Stone Extractor for ureteral stone extraction.
- › The design allows you to engage, reposition, release and extract stones of various sizes throughout the urinary tract.
- › The loop of the basket is highly flexible, kink resistant and durable.
- › A combination that makes it the optimal stone retrieval device for multiple clinical application.

STONE BASKET

Used for manipulation and removal / retrieval of ureteral stones from the bladder, renal pelvis & ureter.

Polyamide sheath basket also available on request:



Sizes

In Fr: Length (in cms):



3 / 4 / 5

70 / 90 / 110

Key Features:

- › Four and five wire configuration are available.
- › Adjustable knob allows for controlled and assured basket expansion.
- › Handle permit complete disassembles.
- › Nitinol baskets have excellent memory & strength.
- › The flat wire construction provides unequaled strength & security.
- › The flat wire basket is designed to more securely grasp and manipulate calculi in the ureter as well as the renal calyx.

PERK BASKET



Sizes

In Fr:	Length (in cms):
8 / 10 / 12	38 / 40 / 42

ZERO TIP

Zero tip is designed for easier maneuverability in the kidney and the ureter.

Sizes

In Fr:	Length (in cms):
1.9 / 2.2 / 3 / 4 / 5	70 / 90 / 110



URS FORCEPS & PCNL FORCEPS



URS Forceps - Rigid & Flexible

Sizes

In Fr:	Length (in cms):
3 / 4 / 5	62

PCNL Forceps - Rigid & Flexible

Sizes

In Fr:	Length (in cms):
7 / 9 / 12	40

Flexible Tri-prong Grasper

Sizes

In Fr:	Length (in cms):
3 / 4 / 5	70 / 90 / 110

Key Features:

- › Tri prong atraumatic forceps
- › Maximum grasping efficiency.

• Also available : Alligator, Peanut Grasper & Stent Removal

NEEDLE

TRU CUT BIOPSY - SEMI AUTOMATIC



Applications

Tru-cut Biopsy needle: A disposable needle with outer cannula and inner notched rod in which a tissue specimen is cut, trapped and withdrawn. The needle is commonly used for percutaneous biopsy of lymph nodes, liver, kidney etc.

Sizes

In Gauge	14G	16G	18G
Length (in cms):	10/ 15 / 20 / 25		

IP NEEDLE



Applications

Used for Percutaneous Nephrostomy puncture. Deep tissue access and all other Percutaneous application.
Used for Nephrostomy Percutaneous application.

Sizes

In Gauge	18G	20G	22G
Length (in cms):	15/ 20/ 23		

Available in : Two part, diamond tip (20G, 22G) / Three part, bevel tip, diamond tip (18G)

CHIBA NEEDLE



Applications

Can be used as a universal needle for percutaneous nephrostomy procedures, cytological biopsies, needle-aspirations, thoracentesis exploratory punctures, to inject contrast dyes etc.

Sizes

In Gauge	18G	20G	22G	23G
Length (in cms):	15, 20, 22 cm			

TURP ELECTRODES 24 FR

Sizes

TURP electrodes compatible with karl storz, Richard Wolf, Olympus working

DOUBLE STEM CUTTING LOOP

DOUBLE STEM ROLLER (3 / 5mm)

SINGLE STEM ROLLER (3 / 5mm)

DOUBLE STEM KOLLINS KNIFE

SINGLE STEM KOLLINS KNIFE

DOUBLE STEM COLD KNIFE



PATH FINDER



60 cc

ASEPTO BULB SYRINGE



12F - 30F / 17cm - 20cm Length

AMPLATZ SHEATH



BLADDER EVACUATOR



60 cc

TOOMY SYRINGE



MALE URINARY PENILE INCONTINENCE CLAMPS



UMBILICAL CLAMP



ALLYX BOTTLE




5X9 mm - 5 meters
8X12 mm - 5 meters

SILICONE TUBING

BILIARY STENT

BILIARY STENT



Applications

Biliary stents provide bile drainage from the gallbladder, pancreas and bile ducts to the duodenum in conditions such as ascending cholangitis due to obstructing gallstones.


Types: Double Pigtail, Single Pigtail, Amsterdam, Straight, Standard, Pancreatic

Sizes

Size (in Fr):	7 FR / 10 FR
Length (in cms):	5 / 7 / 9 / 10 / 12 / 15

CATHETER

PTBD CATHETER SET



Applications

Percutaneous Transhepatic Biliary Drainage (PTBD) Catheter & Set Straight Angled with Pigtail-with or without hydrophilic coating is used for percutaneous drainage applications of biliary and abscess.

Sizes

Size (in Fr):	6 FR - 18 FR
Length (in cms):	30 cm - 60 cm

IUI CATHETER



Applications

The basis behind the IUI procedure is to use a catheter to deliver washed spermatozoa past the cervical mucus barrier into the uterine cavity and thereby increase the sperm concentration at the site of fertilization.

Types: Straight, Angled, Open end & Close end.

Sizes

Size (in Fr):	8FR / 6FR
Length (in cms):	10cm - 20cm

EMBRYO TRANSFER CATHETER




Applications

Embryo transfer is a quick and simple procedure with no discomfort to the patient but also one of the most important steps in the IVF Lifecycle.

Sizes

Size (in Fr):	5FR / 7FR
Length (in cms):	23cm - 25cm

HEMODIALYSIS CATHETER




Applications

Used for exchanging blood to and from the Hemodialysis machine from the patient.

Sizes

Size (in Fr):	9 / 11.5 / 12
Length (in cms):	12cm - 30cm

ENDOTRACHEAL TUBE




Applications

Endotracheal tube is a tube constructed of polyvinyl chloride (PVC) that is placed between the vocal cords through the trachea to provide oxygen and inhaled gases to the lungs.

Sizes

Size (in Fr):	3FR - 32FR
Length (in cms):	20cm / 25cm / 30cm

NBDC CATHETER



Applications

It is used for Biliary duct drainage through nasal passages also for treatment of strictures in the ductile system. It includes a nasal transfer tube and a connecting tube which facilitate easy & safe irrigation and drainage.

Sizes

Size (in Fr):	5 / 6 / 7 / 8 / 9 / 10 / 12
Length (in cms):	180cm - 250cm

NEEDLES

AUTOMATIC BIOPSY NEEDLE



Applications

Automatic biopsy needle: A disposable needle with outer cannula and inner notched rod in which a tissue specimen is cut automatically, trapped and withdrawn. The needle is commonly used for percutaneous biopsy of lymph nodes, liver, kidney, breast & lungs etc.

Sizes

In Gauge	14G	16G	18G
Length (in cms):	10 / 15 / 20 / 25		

SPINAL NEEDLE



Applications

Spinal needles are used to inject analgesia and/or an aesthetic directly into the CSF usually at a point below the second lumbar vertebra. Spinal needles enter the cerebral spinal fluid (CSF) through the membranes surrounding the spinal cord.

Sizes

In Gauge	16G - 24G
Length (in cms):	9cm - 15cm

SCLEROTHERAPY NEEDLE



Applications

Sclerotherapy needles are used for injection of sclerotherapeutic solution in to esophageal varices. Also used for injection of appropriate vasoconstriction/medication agent into the gastrointestinal tract

Sizes

In Gauge	21G / 22G / 23G / 24G / 25G
Length (in cms):	80cm - 250cm

BONE MARROW NEEDLE



Applications

These devices are manual, sterile, disposable needles intended to obtain bone marrow aspirate and core biopsy samples from bone and/or bone marrow.

Sizes

In Gauge	14G / 16G / 18G
Length (in cms):	5cm - 15cm

ENDOMETRIAL TUBE BIOPSY NEEDLE




Applications

Endometrial Tube is used to obtain histological tissues biopsy & extract sample of uterine menstrual content (Sample of Endometrial Tissues).

Sizes

Size (in Fr):	10 FR
Length (in cms):	20cm

OVUM PICKUP NEEDLE



Applications

Used for laparoscopic or ultrasound guided trans vaginal aspiration and flushing of oocytes from ovarian follicles.

Sizes

In Gauge	16G - 25G
Length (in cms):	60cm - 200cm



CONTACT US



SPG



SPG MEDISAFE PRIVATE LIMITED
Imagination | Innovation | Research

- ✉ info@spgmedisafe.com
- ✉ care@spgmedisafe.com
- 🌐 www.spgmedisafe.com



- 📍 **Corporate Address :**
A239, Balaji Bhavan,
Plot No. 42A, Sec - 11, CBD Belapur,
Navi Mumbai- 400614.
Maharashtra, India.
- 📍 **Unit Address :**
SPG Medisafe Pvt. Ltd.,
Platinum Springs,
Unit - No. B-76,77 &78, Plot -6, MIDC,
Taloja- 410208, Maharashtra, India.
- 📞 **Customer Care :**
+91 70455 96408 / 70455 96409

Follow us

-  [spgmedisafe](https://www.linkedin.com/company/spgmedisafe)
-  [spgmedisafe](https://twitter.com/spgmedisafe)
-  [spgmedisafe](https://www.instagram.com/spgmedisafe)
-  [spgmedisafe](https://www.facebook.com/spgmedisafe)

UROSAFE™



Efficient • Transparent • Inclusive